

RAB38 antibody



Catalog Number: 114436

Product name

RAB38 antibody

Immunogen

[Human RAB38 Recombinant protein \(GST tag\)](#)

Specificity

human, mouse; other species not tested.

Antibody description

RAB38 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 27 kDa

Preparation

This antibody was obtained by immunization of RAB38 recombinant protein (Accession Number: NM_022337). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114436(RAB38 antibody) at dilution of 1:500