

Anti-LILRA5 antibody

Catalog Number: 176436

Product name

Anti-LILRA5 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to LILRA5

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human lilra5. aa range:141-190

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of HeLa cells using CD85f Polyclonal Antibody.. Secondary antibody [catalog#: HA1001) was diluted at 1:20000