

Anti-CKM (6A12) Mouse antibody



Catalog Number: 168436

Product name

Anti-CKM (6A12) Mouse antibody

Specificity

Human

Antibody description

CKM (6A12) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human CK-MM protein expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG2b

Applications

ELISA

Validations



Indirect ELISA assay for anti-CK-MM mouse mAb. CKBB as control antigen. Antigen coating concentration: 4ug/ml.