

Sus scrofa (Pig) IL-6 (pig) recombinant protein



Catalog Number: 500836

General Information

Gene Name Synonym

Interleukin-6; Interleukin-6 protein

Protein Construction

A DNA sequence encoding the sus scrofa IL6 (NP_999564.1) (Pro29-Met212) was expressed with an initial Met.

Organism

Sus scrofa (Pig)

Expression Host

E. coli

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant sus scrofa IL6 consists 185 amino acids and predicts a molecular mass of 21.2 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C.

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Sus scrofa (pig) IL6 / IL-6 Protein