

Anti-NT5C1A antibody



Catalog Number: 176136

Product name

Anti-NT5C1A antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to NT5C1A

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human nt5c1a. aa range:151-200

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western blot analysis of lysates from MCF-7 and HepG2 cells, using NT5C1A Antibody. The lane on the right is blocked with the synthesized peptide.




Fig2.; Immunofluorescence analysis of A549 cells, using NT5C1A Antibody. The picture on the right is blocked with the synthesized peptide.