

TIMM9 antibody

Catalog Number: 116136

Product name

TIMM9 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TIMM9 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in A549 cells, mouse liver tissue. Observed molecular weight by Western-blot: 10 kDa

Preparation

This antibody was obtained by immunization of TIMM9 recombinant protein (Accession Number: NM_001304485). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:116136(TIMM9 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:116136(TIMM9 antibody) at dilution of 1:100 (under 10x lens)



Immunofluorescent analysis of HepG2 cells, using

TIMM9 antibody



Catalog Number: 116136

TIMM9 antibody Catalog No:116136 at 1:25

dilution and Rhodamine-labeled goat anti-rabbit
IgG (red).