

Anti-VEGR Receptor 1 Rabbit antibody



Catalog Number: 167136

Product name

Anti-VEGR Receptor 1 Rabbit antibody

Specificity

Human

Antibody description

VEGR Receptor 1 Rabbit polyclonal antibody

Preparation

Antigen: Peptide sequence around aa.31~35 (D-P-E-L-S) derived from Human VEGFR1.

Formulation

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB:1/1000

Validations



Western blot analysis of extract from HeLa, MCF7 and PC12 cells using VEGFR1 rabbit pAb (dilution 1:1000). Predicted band size: 180kDa. Observed band size: 180kDa.