

Anti-MEF2A Mouse antibody



Catalog Number: 169136

Product name

Anti-MEF2A Mouse antibody

Specificity

Human

Antibody description

MEF2A Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human MEF2A (AA: 391-497) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide and 0.5% protein stabilizer

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

ICC/IF, FC

Validations




Western blot analysis using MEF2A mAb against



human MEF2A (AA

Immunofluorescence analysis of HepG2 cells using MEF2A mouse mAb (green). Blue



Flow cytometric analysis of HepG2 cells using MEF2A mouse mAb (green) and negative control

Anti-MEF2A Mouse antibody



Catalog Number: 169136

Western blot analysis using MEF2A mAb against HEK293 (1) and MEF2A (AA

