

CD7 antibody

Catalog Number: 109136

Product name

CD7 antibody

Specificity

Human; other species not tested.

Antibody description

CD7 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human tonsil tissue.

Preparation

This antibody was obtained by immunization of CD7 recombinant protein (Accession Number: NM_006137). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109136(CD7 Antibody) at dilution of 1:200 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109136(CD7 Antibody) at dilution of 1:200 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).