

Anti-CDC42EP1 antibody



Catalog Number: 176536

Product name

Anti-CDC42EP1 antibody

Specificity

Human

Antibody description

Rabbit monoclonal antibody to CDC42EP1

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*TBS (pH7.4), 0.05% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:1,000-1:2,000

Validations




Fig1: Western blot analysis of CDC42EP1 on different lysates using anti-CDC42EP1 antibody at 1/1,000 dilution.; Positive control.; Lane 1: HUVEC; Lane 2: HeLa