

# UBE2S antibody



Catalog Number: 116536

## Product name

UBE2S antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

UBE2S Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, COLO 320 cells, MCF7 cells, MDA-MB-453s cells, mouse ovary tissue. Observed molecular weight by Western-blot: 28kd

## Preparation

This antibody was obtained by immunization of UBE2S recombinant protein (Accession Number: NM\_014501). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:116536(UBE2S antibody) at dilution of 1:800