

Mouse St6galnac2/Siat7/Siat7b (His tag) recombinant protein



Catalog Number: 518536

General Information

Gene Name Synonym

Gal-beta-1,3-GalNAc alpha-2,6-sialyltransferase, GalNAc alpha-2,6-sialyltransferase II, ST6GalNAc II, Sialyltransferase 7B, SIAT7-B

Protein Construction

A DNA sequence encoding the Mouse St6galnac2/Siat7/Siat7b Ser28-Arg373 is expressed with a 6His tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 50mM Tris, 150mM NaCl, pH7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Mouse St6galnac2/Siat7/Siat7b (His tag)
recombinant protein