

CRABP1 antibody



Catalog Number: 109536

Product name

CRABP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CRABP1 Rabbit Polyclonal antibody. Positive WB detected in human spleen tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 15kd

Preparation

This antibody was obtained by immunization of CRABP1 recombinant protein (Accession Number: NM_004378). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:300-1:1500

Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:109536(CRABP1 antibody) at dilution of 1:500