

# Human PPCDC/COAC/MDS018/UNQ9365/PRO34154 (His tag) recombinant protein



Catalog Number: 519536

## General Information

### Gene Name Synonym

Phosphopantothenoylcysteine Decarboxylase,  
PPC-DC, PPCDC, COAC

### Protein Construction

A DNA sequence encoding the Human  
PPCDC/COAC/MDS018/UNQ9365/PRO34154  
Met1-Ser204 is expressed with a 6His tag at the  
N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL  
method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2 µm filtered solution of 20mM  
TrisHCl, 50mM NaCl, 1mM DTT, 10% Glycerol,  
pH 8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human  
PPCDC/COAC/MDS018/UNQ9365/PRO34154 (His  
tag) recombinant protein