

PLAGL1 antibody



Catalog Number: 113936

Product name

PLAGL1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PLAGL1 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, mouse heart tissue. Observed molecular weight by Western-blot: 51 kDa

Preparation

This antibody was obtained by immunization of PLAGL1 recombinant protein (Accession Number: NM_001080953). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:113936(PLAGL1 Antibody) at dilution of 1:600