

Anti-Acetyl-p53 (Lys370) Rabbit antibody



Catalog Number: 166936

Product name

Anti-Acetyl-p53 (Lys370) Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Acetyl-p53 (Lys370) Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human p53 (acetyl K370)

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/500-1/2000

ICC/IF: 1/50-1/200

IP: 1/50

Validations



Western blot analysis of p53 (acetyl K370) expression in HeLa cell lysate treated with Trichostatin A.