

# EME1 antibody



Catalog Number: 110236

## Product name

EME1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

EME1 Rabbit Polyclonal antibody. Positive WB detected in A375 cells, human liver tissue. Observed molecular weight by Western-blot: 65 kDa

## Preparation

This antibody was obtained by immunization of EME1 recombinant protein (Accession Number: NM\_001166131). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:110236(EME1 antibody) at dilution of 1:500