

# DBC1 antibody



Catalog Number: 107236

## Product name

DBC1 antibody

## Specificity

Human, Pig, Mouse; other species not tested.

## Antibody description

DBC1 Mouse Monoclonal antibody. Positive WB detected in fetal human brain tissue, mouse brain tissue, pig brain tissue. Observed molecular weight by Western-blot: 88 kDa

## Preparation

This antibody was obtained by immunization of DBC1 recombinant protein (Accession Number: BC021560). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2c

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107236(DBC1 Antibody) at dilution of 1:1000