

# Human NCR3/1C7/LY117 (Fc tag) recombinant protein



Catalog Number: 518236

## General Information

### Gene Name Synonym

Natural Cytotoxicity Triggering Receptor 3, Activating Natural Killer Receptor p30, Natural Killer Cell p30-Related Protein, NK-p30, NKp30, CD337, NCR3, 1C7, LY117

### Protein Construction

A DNA sequence encoding the Human NCR3/1C7/LY117 Leu19-Thr138 is expressed with a Fc tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 12 months at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  as lyophilized powder. Aliquots of reconstituted samples are stable at  $< -20^{\circ}\text{C}$  for 3 months.

## Formulation

Lyophilized from a  $0.2\ \mu\text{m}$  filtered solution of  $20\text{mM PB}$   $\square$   $150\text{mM NaCl}$ ,  $\text{pH}7.4$ .

## Usage Guide


### Storage

Lyophilized protein should be stored at  $< -20^{\circ}\text{C}$ , though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at  $4-7^{\circ}\text{C}$  for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than  $100\ \mu\text{g}/\text{ml}$ . Dissolve the lyophilized protein in  $\text{ddH}_2\text{O}$ . Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Human NCR3/1C7/LY117 (Fc tag) recombinant protein