

Human PRKCE/PKCE (His tag) recombinant protein



Catalog Number: 519236

General Information

Gene Name Synonym

Protein Kinase C Epsilon Type, nPKC-Epsilon, PRKCE, PKCE

Protein Construction

A DNA sequence encoding the Human PRKCE/PKCE Gln580-Pro737 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PRKCE/PKCE (His tag) recombinant protein