

STAC antibody



Catalog Number: 115636

Product name

STAC antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

STAC Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, human liver tissue. Observed molecular weight by Western-blot: 45kd

Preparation

This antibody was obtained by immunization of STAC recombinant protein (Accession Number: NM_003149). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115636(STAC antibody) at dilution of 1:500