

C11orf42 antibody



Catalog Number: 108636

Product name

C11orf42 antibody

Specificity

Human; other species not tested.

Antibody description

C11orf42 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 36 kDa

Preparation

This antibody was obtained by immunization of C11orf42 recombinant protein (Accession Number: NM_173525). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IF: 1:50-1:500

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108636(C11orf42 Antibody) at dilution of 1:600