

# PUF60 antibody



Catalog Number: 114336

## Product name

PUF60 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PUF60 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, mouse testis tissue, rat testis tissue. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 60 kDa

## Preparation

This antibody was obtained by immunization of PUF60 recombinant protein (Accession Number: NM\_014281). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114336(PUF60 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:114336(PUF60 antibody) at dilution of 1:100 (under 10x lens)