

# Transketolase antibody



Catalog Number: 116336

## Product name

Transketolase antibody

Recommended Dilution:

WB: 1:200-1:2000

## Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

IF: 1:10-1:100

## Antibody description

Transketolase Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, BxPC-3 cells, HEK-293 cells, human kidney tissue, Jurkat cells, MCF7 cells, mouse brain tissue. Positive IHC detected in human lymphoma tissue, human liver tissue. Positive IF detected in BxPC-3 cells. Observed molecular weight by Western-blot: 68 kDa

## Validations




## Preparation

This antibody was obtained by immunization of Transketolase recombinant protein (Accession Number: NM\_001064). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116336(TKT antibody) at dilution of 1:300

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:116336(TKT antibody) at dilution of 1:100 (under 10x lens)

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IF


## Dilutions

# Transketolase antibody



Catalog Number: 116336

---



Immunofluorescent analysis of BxPC-3 cells using Catalog No:116336(TKT Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)