

Anti-CD166/ALCAM antibody



Catalog Number: 100736

Product name

Anti-CD166/ALCAM antibody

Immunogen

[Cynomolgus CD166/ALCAM \(His Tag\) recombinant protein](#)

Specificity

Cynomolgus CD166 / ALCAM

Antibody description

Rabbit monoclonal to CD166/ALCAM

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Cynomolgus CD166 / ALCAM (G7NZQ8; Met1-Ala526).

Formulation

0.2 µm filtered solution in PBS with 5% treha

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 5-10 µg/ml

Validations



Items	Lanes	A	B
Sample (whole cell lysate)		HuT78	Daudi
Sample Volume (µg/lane)		30	30
Gel	13% SDS-PAGE reducing gel		
Recommended Concentration	5-10 µg/ml		
Secondary Antibody	Dylight 800 labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.		
Developed using Odyssey imaging system.			
Explanation	Predicted band size : 65 kDa Observed band size : 65 kDa		

CD166 / ALCAM Antibody, Rabbit MAb, Western blot