

SLC12A6 Polyclonal Antibody



Catalog Number: 165736

Product name

SLC12A6 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to SLC12A6

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 60-190 of human SLC12A6 (NP_598408.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - SLC12A6 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using SLC12A6 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.