

# Anti-BRCC36 Rabbit antibody



Catalog Number: 166736

## Product name

Anti-BRCC36 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

BRCC36 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human BRCC36

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, FC, IP


## Dilutions

WB: 1/2000-1/10000

FC: 1/20-1/100

IP: 1/20-1/50

## Validations



Western blot detection of BRCC36 in Rat Brain,C6,3T3,Hela cell lysates using BRCC36 Rabbit pAb(1:1000 diluted).Predicted band size:36kDa.Observed band size:36kDa.