

Catalog Number: 108736

Product name

C5orf26 antibody

Specificity

Human; other species not tested.

Antibody description

C5orf26 Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer tissue, human brain tissue. Positive IF detected in HepG2 cells, HeLa cells.

Preparation

This antibody was obtained by immunization of C5orf26 recombinant protein (Accession Number: BC112077). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:10-1:100


Validations



Immunohistochemistry of paraffin-embedded human liver cancer slide using Catalog No:108736(C5orf26 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human liver cancer slide using Catalog No:108736(C5orf26 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:108736(C5orf26 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG