

# Anti-Cytokeratin 9 Rabbit antibody



Catalog Number: 168046

## Product name

Anti-Cytokeratin 9 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Cytokeratin 9 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human Cytokeratin 9

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, IP


## Dilutions

WB: 1/1000

IHC: 1/50

IP: 1/20

## Validations



Western blot detection of Cytokeratin 9 in C6,Ramos using Cytokeratin 9 Rabbit mAb(1:1000 diluted)