

Anti-Factor D Rabbit antibody



Catalog Number: 169046

Product name

Anti-Factor D Rabbit antibody

Specificity

Human

Antibody description

Factor D Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Factor D

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

1:1000

Validations



Western blot detection of Factor D in HeLa cell lysates using Factor D Rabbit mAb(1:1000 diluted).Predicted band size:27kDa.Observed band size:27kDa.