

# Anti-Factor D Rabbit antibody



Catalog Number: 169046

## Product name

Anti-Factor D Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB

## Antibody description

Factor D Rabbit monoclonal antibody

## Dilutions

1:1000

## Preparation

Antigen: Recombinant protein of human Factor D

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Validations



Western blot detection of Factor D in HeLa cell lysates using Factor D Rabbit mAb(1:1000 diluted). Predicted band size: 27kDa. Observed band size: 27kDa.