

TSPAN15 antibody



Catalog Number: 116446

Product name

TSPAN15 antibody

Specificity

Human; other species not tested.

Antibody description

TSPAN15 Rabbit Polyclonal antibody. Positive WB detected in human ovary cancer. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 33kd

Preparation

This antibody was obtained by immunization of TSPAN15 recombinant protein (Accession Number: NM_012339). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human ovary cancer were subjected to SDS PAGE followed by western blot with Catalog No:116446(TSPAN15 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:116446(TSPAN15 antibody) at dilution of 1:100 (under 10x lens)