

# TALPID3 antibody



Catalog Number: 115846

## Product name

TALPID3 antibody

## Immunogen

[Human KIAA0586 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TALPID3 Rabbit Polyclonal antibody. Positive WB detected in mouse embryo tissue. Observed molecular weight by Western-blot: 60 kDa

## Preparation

This antibody was obtained by immunization of TALPID3 recombinant protein (Accession Number: NM\_014749). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



mouse embryo tissue were subjected to SDS PAGE followed by western blot with Catalog No:115846(KIAA0586 Antibody) at dilution of 1:600