

CAPZB antibody

Catalog Number: 108846

Product name

CAPZB antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CAPZB Rabbit Polyclonal antibody. Positive IF detected in A431 cells. Positive WB detected in Jurkat cells, A431 cells, HEK-293 cells, HeLa cells, mouse testis tissue. Observed molecular weight by Western-blot: 30-31 kDa

Preparation

This antibody was obtained by immunization of CAPZB recombinant protein (Accession Number: NM_004930). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108846(CAPZB Antibody) at dilution of 1:1000



Immunofluorescent analysis of (-20oc Ethanol) fixed A431 cells using Catalog No:108846(CAPZB Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)