

GPT antibody

Catalog Number: 111146

Product name

GPT antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GPT Rabbit Polyclonal antibody. Positive IP detected in L02 cells. Positive WB detected in human heart tissue, mouse kidney tissue, mouse liver tissue, rat liver tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 50-55 kDa

Preparation

This antibody was obtained by immunization of GPT recombinant protein (Accession Number: NM_005309). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111146(GPT antibody) at dilution of 1:400



IP Result of anti-GPT (IP:Catalog No:111146, 4ug; Detection:Catalog No:111146 1:500) with L02 cells lysate 1800ug.