

# ARL8A antibody

Catalog Number: 108146

## Product name

ARL8A antibody

## Immunogen

[Human ARL8A Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ARL8A Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, HepG2 cells. Observed molecular weight by Western-blot: 22 kDa

## Preparation

This antibody was obtained by immunization of ARL8A recombinant protein (Accession Number: NM\_138795). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108146(ARL8A Antibody) at dilution of 1:600



IP Result of anti-ARL8A (IP:Catalog No:108146, 4ug; Detection:Catalog No:108146 1:600) with mouse brain tissue lysate 3440ug.