

CD8b antibody

Catalog Number: 109146

Product name

CD8b antibody

Specificity

Human; other species not tested.

Antibody description

CD8b Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 37-43 kDa

Preparation

This antibody was obtained by immunization of CD8b recombinant protein (Accession Number: NM_172213). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109146(CD8B antibody) at dilution of 1:300