

RAP1GAP antibody



Catalog Number: 114546

Product name

RAP1GAP antibody

Immunogen

[Human RAP1GAP Recombinant protein \(His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

RAP1GAP Rabbit Polyclonal antibody. Positive IF detected in HeLa. Positive WB detected in HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 85-95kd

Preparation

This antibody was obtained by immunization of RAP1GAP recombinant protein (Accession Number: XM_017001988). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114546(RAP1GAP antibody) at dilution of



Immunofluorescent analysis of HeLa cells, using RAP1GAP antibody Catalog No:114546 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).