

# Anti-DOHH Rabbit antibody



Catalog Number: 167546

## Product name

Anti-DOHH Rabbit antibody

## Specificity

Human

## Antibody description

DOHH Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human DOHH

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions

WB: 1/1000

IP: 1/20

## Validations



Western blot detection of DOHH in K562, HeLa cell lysates using DOHH Rabbit pAb(1:1000 diluted). Predicted band size: 33kDa. Observed band size: 33kDa.