

BTN2A2 antibody



Catalog Number: 108546

Product name

BTN2A2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

BTN2A2 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, NIH/3T3 cells, human testis tissue. Observed molecular weight by Western-blot: 59kd

Preparation

This antibody was obtained by immunization of BTN2A2 recombinant protein (Accession Number: XM_005248797). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108546(BTN2A2 antibody) at dilution of 1:500