

Human SULT2A1/HST/STD (His tag) recombinant protein



Catalog Number: 518546

General Information

Gene Name Synonym

Bile Salt Sulfotransferase,
Dehydroepiandrosterone Sulfotransferase, DHEA-ST, Hydroxysteroid Sulfotransferase, HST, ST2, ST2A3, Sulfotransferase 2A1, ST2A1, SULT2A1, HST, STD

Protein Construction

A DNA sequence encoding the Human SULT2A1/HST/STD Ser2-Glu285 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 100mM NaCl, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human SULT2A1/HST/STD (His tag) recombinant protein