

GFI1 antibody



Catalog Number: 110946

Product name

GFI1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GFI1 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HL-60 cells, PC-3 cells, Raji cells, RAW264.7 cells. Observed molecular weight by Western-blot: 62-66 kDa

Preparation

This antibody was obtained by immunization of GFI1 recombinant protein (Accession Number: NM_001127216). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110946(GFI1 antibody) at dilution of 1:500