

PLCD1 antibody



Catalog Number: 113946

Product name

PLCD1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PLCD1 Rabbit Polyclonal antibody. Positive WB detected in human rectum tissue, human brain tissue. Observed molecular weight by Western-blot: 85kd

Preparation

This antibody was obtained by immunization of PLCD1 recombinant protein (Accession Number: NM_006225). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human rectum tissue were subjected to SDS PAGE followed by western blot with Catalog No:113946(PLCD1 antibody) at dilution of 1:400