

# AKR7A3 antibody



Catalog Number: 107946

## Product name

AKR7A3 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

AKR7A3 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, HepG2 cells, human brain tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 37kd and 55-60kd

## Preparation

This antibody was obtained by immunization of AKR7A3 recombinant protein (Accession Number: NM\_012067). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107946(AKR7A3 antibody) at dilution of 1:400