

Anti-PAX1 antibody



Catalog Number: 176246

Product name

Anti-PAX1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to PAX1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human pax-1. aa range:318-367

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of 293T cells using Pax-1 Polyclonal Antibody diluted at 1: 500 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).