

# PECI antibody

Catalog Number: 113646

## Product name

PECI antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

PECI Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, human heart tissue, human placenta tissue, Jurkat cells. Positive IHC detected in human heart. Observed molecular weight by Western-blot: 40kd

## Preparation

This antibody was obtained by immunization of PECI recombinant protein (Accession Number: NM\_001166010). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113646(PECI antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human heart using Catalog No:113646(PECI antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:113646(PECI antibody) at dilution of 1:100 (under 40x lens)