

Sprr2i antibody



Catalog Number: 115646

Product name

Sprr2i antibody

Specificity

Human; other species not tested.

Antibody description

Sprr2i Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue.

Preparation

This antibody was obtained by immunization of Sprr2i recombinant protein (Accession Number: NM_011475). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:115646(Sprr2i Antibody) at dilution of 1:50 (under 40x lens)