

# Human BCL2L1/BCL2L/BCLX (His tag) recombinant protein



Catalog Number: 518346

## General Information

### Gene Name Synonym

Bcl-2-Like Protein 1, Bcl2-L-1, Apoptosis Regulator  
Bcl-X, BCL2L1, BCL2L, BCLX

### Protein Construction

A DNA sequence encoding the Human  
BCL2L1/BCL2L/BCLX Met1-Arg212 is expressed  
with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL  
method.

## Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

## Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM  
HEPES, 50mM KCl, 20% Glycerol, pH 7.5.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

## SDS-PAGE



Human BCL2L1/BCL2L/BCLX (His tag)  
recombinant protein