

CLEC10A/CD301 antibody



Catalog Number: 109346

Product name

CLEC10A/CD301 antibody

Specificity

Human; other species not tested.

Antibody description

CLEC10A/CD301 Rabbit Polyclonal antibody.
Positive WB detected in BxPC-3 cells. Observed
molecular weight by Western-blot: 32 kDa, 28 kDa

Preparation

This antibody was obtained by immunization of
CLEC10A/CD301 recombinant protein (Accession
Number: NM_182906). Purification method:
Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



BxPC-3 cells were subjected to SDS PAGE
followed by western blot with Catalog
No:109346(CLEC10A antibody) at dilution of
1:400