

# H5N1 HA (His Tag) recombinant protein



Catalog Number: 503746

## General Information

### Protein Construction

A DNA sequence encoding the Influenza A virus (A/chicken/VietNam/NCVD-016/2008(H5N1)) hemagglutinin (ACO07033.1) (Met1-Arg343), termed as HA1, was expressed with a polyhistidine tag at the C-terminus.

### Organism

H5N1

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Asp 17

### Molecular Mass

The recombinant HA1 subunit of Influenza A virus (A/chicken/VietNam/NCVD-016/2008(H5N1)) consists 338 amino acids and predicts a molecular

mass of 38.2 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage


Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Influenza A H5N1  
(A/chicken/VietNam/NCVD-016/2008)  
Hemagglutinin Protein (HA1 Subunit) (His Tag)  
SDS-PAGE